



AQUA SONIC

Aqua Therapy System

AQUA SONIC is an advanced aqua therapy device designed to improve skin exfoliation by removing aged dead skin cells and excess sebum while delivering deep hydration. Its uniquely designed tip features two suction holes and one serum discharge hole, with a spiral structure that ensures thorough skin contact for effective softening of dead skin cells. The built-in vacuum function enhances both the peeling and extraction processes.

Additionally, AQUA SONIC provides a range of specialized handpieces to help achieve optimal clinical results.



Diamond Handpiece



Hydrofacial Handpiece



Ultrasound Handpiece



Feature

Aqua therapy

Eliminates aged dead skin cells and excess sebum, and provides hydration at the same time.

Special tip design

Two suction holes and one serum discharge hole, incorporating a spiral design to ensure sufficient contact to skin.

Vacuum function

Amplifies the peeling and extrusion effects.

Curated solutions

Enhances efficacy when paired with specially curated solutions like LASER DOCTOR and other compatible options.

Technical specification

| | |
|-----------------------|---|
| Ultrasound Handpiece | 1.2 MHz |
| Diamond Handpiece | #140 / #170 |
| Hydrofacial Handpiece | Face (regular size), nose (narrow size) |
| Vacuum Pump | Up to 600 mmHg (±10 %) |
| Power Consumption | 80 W |
| Aqua Waste Tank | 1200 cc |
| Input Voltage | 200 VAC / 50 / 60 Hz |
| Dimension (W*L*H) | 350 mm (W) X 350 mm (L) X 1000 mm (H) |
| Weight | 30 Kg |

